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Please read the following instructions carefully before use.

Product specification:

Prodox urine test is a qualitative and semi-quantitative test for urine analysis. It is intended for marking of glucose, bilirubin, ketones, specific gravity, blood, pH, protein, urobilinogen, nitrites and leukocytes in dog's urine.

Contents of the kit:

- 2x pouch with the test strip - individually wrapped
- 1x color card of the results
- 2x tube
- 2x pipette
- 2x single-use gloves
- 1x instruction of use

Collecting the samples:

- Test the sample as soon as possible after collection.
- If the urine test can't be performed within 1h after urine collection, don't perform the test.
- Collect urine coming from the middle stream, preferably from the first morning collection.
- Do not collect the urine faster than 2h after the last meal, too short an interval from a meal will falsify the results.

Interpretation of test results:

Results are obtained by direct comparison of the test strip to the color chart included in the kit. Very high concentrations of analytes may cause unusual color reactions. In such cases it is recommended that users interpret results based on color intensity.

GLUCOSE: Glucose belongs to simple sugars and is used by the body as one of its primary sources of energy. Glucose is not detected in the urine of healthy dogs. Filtered glucose is almost entirely absorbed by the renal tubule, and only a very small amount is excreted in the urine.

BILIRUBIN: Bilirubin is an orange-red bile pigment - the final product of the metabolism of red blood cell hemoglobin (the primary role in this process is played by the liver). Normally, no bilirubin is found in the urine, the exception being male dogs, in which trace amounts may appear.

KETONES: Ketone bodies are compounds that are intermediate products of fat metabolism. In a healthy body, their trace amounts are observed in the blood serum. If there is an increased production of ketone bodies, it is possible to detect these compounds in the urine. Normally, ketone bodies are not found in the urine, their trace amounts can appear as a result of animal's starvation. A positive result occurs with diabetes, fever or liver damage.

SPECIFIC GRAVITY: Specific gravity shows whether the kidneys are properly densifying the urine. Normal values for dogs are 1,016 - 1,045. However the specific gravity should be measured with refractometer. This value should be interpreted in conjunction with the patient's hydration, blood urea nitrogen (BUN) and creatinine levels.

BLOOD: The results of the test for this parameter in urine collected from healthy dogs should be negative. A positive reaction for latent blood indicates the presence of erythrocytes, hemoglobin or myoglobin in the urine. Any positive result should prompt a thorough analysis of the urine sediment via microscopic examination.

pH: Urine pH in animals depends mainly on diet and acid-base balance. Meat and high-protein diets provide an acidic urine reaction (due to excretion of acidic products of protein metabolism), but this can vary in healthy dogs. Normal pH values range between 5,5 - 7,0 in dogs.

PROTEIN: Protein in urine can be qualitatively assessed using urine dip strips. Healthy dogs excrete negligible amounts of protein, so a small level of protein in a randomly gathered sample is not necessarily a cause for concern. In case of a positive result, it is recommended to contact your veterinarian for more detailed tests (such as determining the level of UPC in urine).

UROBILINOGEN: Urobilinogen is a pigment derived from the metabolism of bilirubin. It is formed in the lumen of the gastrointestinal tract, from where approx. 20% is absorbed into the blood and, passing through the liver, excreted in the bile and, in small amounts, also in the urine through the kidneys. Normal urine contains small amounts of urobilinogen. Increased levels may indicate liver damage, digestive tract issues, purulent uteritis or excessive breakdown of erythrocytes. Normal urine contains small amounts of urobilinogen. Increased levels may indicate liver damage, digestive tract issues, purulent uteritis or excessive breakdown of erythrocytes.

NITRITES: Normally nitrites are not detectable in urine. Testing urine for the presence of nitrites may indicate the presence of bacteria in urine. Bacteriuria is usually caused by infection of the kidneys, ureters, bladder, or urethra. However, majority of bacterial infections in dogs is not caused by organisms that reduce nitrate to nitrite. Therefore in veterinary medicine both false positive and false negative results occur, this test has too low sensitivity for widespread use.

LEUKOCYTES: Leukocytes are white blood cells that are involved in immune reactions. An increase in their number accompanies inflammation and infections. Note that leukocytes in urine from micturition can come from both the urinary and reproductive systems. Trace amounts of leukocytes in urine can be found normally. It is not recommended to determine the number of leukocytes in dogs using test strips. False positive or false negative results may occur. Leukocyturia is diagnosed by examining the urine sediment.


Storage conditions:

- Store in room temperature 2°C-30°C, keep away from direct sunlight.
- Keep away from children.
- Do not use after expiration date specified on the label.

Warning and precautions:


- Be careful when performing the test, do all the steps wearing single-use gloves that are part of this kit.
- Do not touch the test area of the strip.
- Protect the test strips from humidity, heat and direct sunlight.
- After use, you can dispose of the contents of the kit in mixed household waste.
- All packaging elements are non-toxic and safe if used as instructed.
- Do not use if the product packaging is damaged.

Prodox urine test strip measures chemical properties of urine, they are only one of the elements of a proper examination. Comprehensive urine examination also includes physical and microscopic testing of urine sediment. As with all diagnostic tests, the final diagnosis should not be based on a single test result but should be established by the veterinarian only after evaluation of all clinical signs and laboratory results.

EN Take the urine sample into a clean, dry, unused container. Put the urine sample into the tube using the pipette. Close the tube with the cap. Remove the Prodox Urine Test strip from the pouch. Take off the cap from the tube.

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EN Immerse all reagent areas in the urine sample for 1-2 seconds. Remove excess urine from the strip by delicately tapping the edge of the strip on the sample container or gently dry the edges of the strip with a tissue. Hold the test strip in a horizontal position. Read the results within 30-60 seconds. Compare the results to the color chart with result interpretation included in the test kit. Color changes that appear only along the edges of the area or after a certain amount of time are not diagnostically significant.

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EN After use collect all the test components to an empty trash bag, close it precisely. Throw it out to mixed trash.

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PRODOX
Professional Vital PetCare

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Prodox Vital GmbH